



140626 Ansel Lab/ Sana Patel

- 30  $\mu\text{L}$  rxn is now at 24ng/ $\mu\text{L}$  (of EACH "set" of phosphorylated oligo [2 sets])
- Dilute phosphorylated oligo mix: 5 $\mu\text{L}$  phosphorylated mix + 15 $\mu\text{L}$  MQ water;  $\rightarrow$  6ng/ $\mu\text{L}$  of each oligo

**e) Not1/EcoRV digest & dephosphorylate pEntr-GFP vector** (If you already have a cut/dephosphorylated pEntrGFP vector, skip to step f)

- 1) Digest pEntr-GFP vector with Not1 and EcoRV. Example rxn (change based on how much DNA you have):
  - 2.5uL EcoRV (20,000U/mL)
  - 5 uL Not I (10,000U/mL)
  - 7.5 uL 10x NEBuffer 3.1
  - X uL pEntr-GFP (5ug)
  - add H2O to 75uL
- 2) Run digestion mix in 1% agarose gel.
- 3) Gel extract the digested band (~2.8kb). (we use Sigma Gel extract kit)
- 4) Dephosphorylate the cut/purified vector with Antarctic Phosphatase (NEB) so it doesn't ligate to itself, 30min-1hr at 37C. Example rxn:
  - 25 uL of cut/purified vector
  - 2uL Antarctic Phosphatase
  - 5uL 10X antarctic phosphotase buffer
  - 18uL H2O
- 5) PCR clean up the dephosphorylated vector (Sigma kit) to get rid of buffer, enzymes and elute the DNA in about 40 $\mu\text{L}$ .
- 6) Nanodrop to measure concentration. Dilute/concentrate cut vector to 50ng/ $\mu\text{L}$

**f) Ligate sensor sequence into pEntr-GFP vector**

- We use Quick Ligation kit from NEB (M2200)

- 1) Combine 50ng of vector with a 3-fold and 6-fold molar excess of insert to a final volume of 10uL. Include a "no insert" control reaction. Concept:
  - pEntr: ~2800bp
  - inserts: ~60bp each = approx 50x smaller than vector
  - Vector's MW is approx 50x as much as the insert's because it's approx 50x as long, so you need 50x more vector than insert for a 1:1 molar ratio (to get the same number of "ends" to ligate). For a 1:1 molar ratio with 50ng vector, you would need 1ng of insert. For a 3:1 ratio, you need 3ng of insert. 6:1 = 6ng of insert.

Example 6-fold molar excess reaction:

- 1 uL pEntr-GFP cut vector (50ng/uL)
  - 1 uL of phosphorylated sensor oligos (6 ng/uL)
  - 8 uL MQ H<sub>2</sub>O
- 2) Add 10 $\mu\text{L}$  of 2x Quick ligation Reaction buffer
  - 3) Add 1 $\mu\text{L}$  of Quick T4 DNA ligase and mix thoroughly. Quick vortex.
  - 4) Centrifuge briefly and incubate at RT for 5min.
  - 5) Chill on ice, then transform bacteria or store at -20C.

**g) Transform ligated pEntr-GFP-sensor plasmid into competent cells.**

- 1) Thaw competent cells on ice. (We use NEB Top10 competent cells).
- 2) Each vial has ~50 $\mu\text{L}$  of bacteria. Split 25uL of bacteria into eppies per ligation reaction, plus one extra for a "No Insert" control.

140626 Ansel Lab/ Sana Patel

- 3) Add 2 $\mu$ L of ligation mix into 25 $\mu$ L of bacteria.
- 4) Incubate on ice for 15-30min.
- 5) Heat shock for 2 min at 37C, or 45 sec on 42C; chill on ice for 5min.
- 6) Add SOC media or LB broth at 9:1, next to flame. (180 $\mu$ L per 25 $\mu$ L of bacteria)
- 7) Incubate in shaker or water bath at 37C for 1hr. (For Kanamycin-resistant, make sure to incubate for at least 1 hr; when doing Ampicillin-selection can cut to 30-45 min.)
- 8) Pre-heat and dry KANAMYCIN-resistant LB plate to 37C, while bacteria is incubating.
- 9) Spread 25-30% (~50 $\mu$ L) of total bacteria mix onto KAN-plate.
- 10) Incubate over-night in 37C bacterial incubator.

**h) Got colonies? Screen them by bacterial PCR.** (If your "no insert" plate has few to no colonies relative to the other plates, skip to h-4, but pick multiple colonies [~3-5] to miniprep and sequence)

Add 100 $\mu$ L of MQ water into 96-well plate (I prefer U-bottom plates with lids). I usually screen at least 12 colonies from each plate. And 3-4 colonies from control "No insert".

1. With a 20 $\mu$ L pipette tip, poke each colony and transfer a small amount of bacteria into each well. Set the pipette at 20 $\mu$ L and pipette up/down several times to flush out the bacteria into the well. Label each colony with the well # it goes to in the 96-well plate.
2. Keep plate on ice. Bacteria will lyse by itself eventually. If not, PCR denaturation step will take care of it. Plate can be stored at -20C.
3. Make a PCR mix:
  - a. 10  $\mu$ L 2x MyTaq Red
  - b. 0.2  $\mu$ L eGFP F (25 $\mu$ M)
  - c. 0.2  $\mu$ L mature miRNA seq (25 $\mu$ M)
  - d. 4.6  $\mu$ L H<sub>2</sub>O

**Primers:** we use Forward primer for GFP and the microRNA mature DNA sequence as reverse primer. The PCR band is ~750bp.

EGFP Fwd primer: AGCTGACCCTGAAGTTCATCTG

- If your pEntr-GFP vector came from digested pEntr-GFP-sensor plasmid, include a separate PCR using the GFP Forward primer + primer for that specific sensor. There might be some undigested sensor in your DNA mix.

- Aliquot 15 $\mu$ L of mix into PCR plate

- With a multi-channel, add 5 $\mu$ L of Bacterial prep into PCR mix.

- e. PCR program:
  - i. 94C 3min
  - ii. 35 cycles: 94C 15sec, 58C 25sec, 72C 45sec
  - iii. 72C 10min
  - iv. Hold at 4C
- f. Load most of PCR mix (maybe 17 $\mu$ L) in 1% agarose gel.
- g. If you get a band for YOUR sensor sequence and not for the other sensor (where pEntr came from)  $\rightarrow$  miniprep and sequence that colony. Do more than 1 if you have multiple positive clones.
4. Grow the right colony, by dipping a 200 $\mu$ L pipette tip into the colony and dropping it into 3mL of LB broth + KANAMYCIN. Grow O/N in shaker at 37C.
5. Miniprep using Sigma kit, have DNA eluted at least at 100ng/ $\mu$ L for sequencing (if less, it's okay, try to sequence anyway).
6. Send for sequencing using Sequetech, select "Plasmid" and use their own cNew eGFP primers.
  - a. When get results back, make sure all 4 microRNA binding sequence are there. Also check it's in right direction and enzymatic sites are correct.
  - b. Got right pEntr-GFP-sensor plasmid? Clonase it into MSCV-hCD25+ccdb vector

**i) Gateway LR clonase reaction:** the whole rxn can be halved (75ng + 75ng in 4 $\mu$ L + 1 $\mu$ L enzyme)

Gateway LR Clonase II: Invitrogen 11791-020

1. Add 50-150ng of entry clone + 150ng MSCV-hCD25+ccdb (150ng). Adjust volume to 8 $\mu$ L with TE buffer, pH8.0. Include a "no entry clone" control reaction.
2. Thaw the LR clonase II enzyme mix on ice for about 2 minutes. Vortex the LR clonase II enzyme, mix briefly twice (2 sec each time).
3. Add 2 $\mu$ L of LR clonase enzyme II to the reaction and mix well by vortexing twice. Quick spin down. Return enzyme to -20C
5. Incubate the clonase reaction at 25C for 1 hr.
6. Add 1 $\mu$ L of Proteinase K solution to terminate reaction. Vortex briefly. Incubate samples at 37C for 10min.
7. Can store in -20C, or transform into bacteria and spread it on AMPICILLIN LB plates.
8. Pick 2 clones  $\rightarrow$  miniprep  $\rightarrow$  diagnostic digest of  $\sim$ 1 $\mu$ g with EcoRV  $\rightarrow$  run on a 1% agarose gel to confirm clonase rxn worked (you should see 3 large bands: >3kb, >2kb, >1kb)

**j) Maxiprep!**

1. Pick the colony confirmed by the EcoRV digest and make a starter culture: drop tip in 3ml LB broth+ AMP. Do this in AM.
2. Incubate in 37C shaker for 8hrs.
3. Add 150-250mL LB broth + ampicillin into Erlenmeyer flask that is 10x larger than the volume of your starter culture
4. Add starter culture into broth at 1:750.
5. Incubate O/N in shaker at 37C.
6. Maxiprep culture following instructions of Sigma Maxiprep kit. Select elution option for maximum recovery of plasmid, and then complete the DNA concentration step.
7. Use a sharpie to mark the location of the DNA pellet, as it will turn clear while the EtOH dries. Resuspend in  $\sim$ 500 $\mu$ L TE. Let sit on ice for 10min to dissolve completely. Pipette over the location of the pellet - if liquid is sticking to the tube there, the DNA has not fully gone into solution.
9. Nanodrop. Dilute to 1 $\mu$ g/ $\mu$ L. Store in -20C.

You can create sensor-bearing virus through Calcium Phosphate Transfection of phoenix cells, and then transduce mouse CD4+ T cells with the resulting virus.